




LIPOSUCTION CANNULAS



Cannulas For Iniltration
Measurements = working length



S26-101-02
2mmØ
25cm
10holes



S26-102-03
3 mm Ø
25 cm
10 holes



Cannulas
S26-103-01
9 x 70 mm Ø
S26-104-01
1 x 120 mm Ø
Packing unit 100 pieces, steril



S26-105-01
Luer-Lock syringes
Adaptor with Luer-Lock connection, for use of
cannulas with handle connector on



S26-106-01
Accessories for Syringes




**Entner
Cannulas**



S26-111-02
2 mm Ø
15 cm
18 holes



S26-112-03
3 mm Ø
15 cm
18 holes



Entner Cannulas
Measurements = working length



S26-113-02
2 mm Ø
15 cm
18 holes



S26-114-03
3 mm Ø
15 cm
18 holes



Entner Cannulas
Measurements = working length



S26-115-03
3 mm Ø
20 cm
22 holes



S26-116-04
4 mm Ø
20 cm
22 holes



**Entner
Cannulas**



S26-117-03
3 mm Ø
20 cm
22 holes



S26-118-04
4 mm Ø
20 cm
22 holes

Liposuction Cannulas
Cánulas de liposucción
Liposuction-Kanülen
Canules pour liposuction




**Entner
Cannulas**


S26-138-03
3 mm Ø
20 cm
22 holes



**Entner
Cannulas**

Measurements = working length

S26-139-03
3 mm Ø
20 cm
22 holes



**Entner
Cannulas**

Measurements = working length

For fine facial and neck
contouring.

S26-140-02
2 mm Ø
15 cm
8 holes



**Entner
Cannulas**

Measurements = working length

S26-141-02
2 mm Ø
15 cm
8 holes

Set Suggestions Liposuction Cannulas
Setvorschläge Liposuction-Kanülen

Proposiciones de juegos cánulas de liposucción
Propositions pour boîtes canules pour liposuction



Face Cannulas
S26-142-02
 15 cm
 2 mm Ø



S26-142-01
 5 cm 1,7 mm Ø
S26-143-01
 15 cm 1,7 mm Ø



Neck Cannulas
S26-144-03
 15 cm
 3 mm Ø



Neck Cannulas
S26-145-03
 20 cm
 3 mm Ø



Neck Cannulas
S26-146-02
 15 cm
 2 mm Ø



S26-147-02
 5 cm 2 mm Ø
S26-148-03
 5 cm 3 mm Ø
S26-149-02
 15 cm 2 mm Ø
S26-150-03
 15 cm 3 mm Ø



Thigh - Lower Leg Cannulas
S26-151-03 20 cm 3 mm Ø
S26-152-03 30 cm 3 mm Ø
S26-153-04 20 cm 4 mm Ø
S26-154-04 30 cm 4 mm Ø



Abdomen / Saddlebag Cannulas
S26-155-03 20 cm 3 mm Ø
S26-156-04 20 cm 4 mm Ø



Abdomen / Saddlebag Cannulas
S26-157-03 30 cm 3 mm Ø
S26-157-04 30 cm 4 mm Ø



S26-158-00
 with rotatable tube connection
 Accessories

Liposuction Cannulas
Cánulas de liposucción
Liposuction-Kanülen
Canules pour liposuction

Fork shape, rounded tips, dissecting surface between tips / one central hole.



S26-161-02
2 mm Ø
5 cm



S26-163-03
3 mm Ø
5 cm



S26-167-04
4 mm Ø
5 cm



S26-170-05
5 mm Ø
30 cm

S26-162-02
2 mm Ø
15 cm

S26-164-03
3 mm Ø
15 cm

S26-168-04
4 mm Ø
25 cm

S26-165-03
3 mm Ø
25 cm

S26-169-04
4 mm Ø
30 cm

S26-166-03
3 mm Ø
30 cm

Toledo Cannulas

Measurements = working length



S26-171-03
3 mm Ø
15 cm



S26-173-04
4 mm Ø
15 cm

S26-172-03
3 mm Ø
25 cm

S26-174-04
4 mm Ø
25 cm

three lateral holes



S26-175-03
3 mm Ø
30 cm

Turbo



The improved arrangement of the suction holes ensures a more effective and gentle suction in areas with much connective tissue and in case of gynecomasty.

Mercedes Cannulas

Measurements = working length

Kotzur Cannulas



Cannulas



one central hole, two lateral holes
Measurements = working length

S26-181-03
3 mm Ø
20 cm

S26-182-03
3 mm Ø
30 cm



S26-183-04
4 mm Ø
30 cm



S26-184-05
5 mm Ø
30 cm



S26-185-06
6 mm Ø
30 cm



Cannulas

Measurements = working length



two lateral holes,
one central hole

S26-186-03
3 mm Ø
20 cm

S26-187-03
3 mm Ø
30 cm



S26-188-04
4 mm Ø
30 cm



S26-189-05
5 mm Ø
30 cm



S26-190-06
6 mm Ø
30 cm





S26-191-03
3 mm Ø
20 cm



S26-192-03
3 mm Ø
30 cm



S26-193-04
4 mm Ø
30 cm



S26-194-04
4 mm Ø
40 cm



S26-195-05
5 mm Ø
30 cm



S26-196-06
6 mm Ø
30 cm



S26-197-04
4 mm Ø
30 cm



S26-198-05
5 mm Ø
30 cm



S26-199-03
3 mm Ø
20 cm



S26-200-03
3 mm Ø
30 cm



S26-201-30
30 cm



Cannulas

three holes, diagonal

Measurements = working length

Cannulas

one central hole, two lateral holes

Measurements = working length

Cannulas

Point with hole, one lateral hole

Cannulas

Measurements = working length

For facial contouring.
one lateral hole



Cannulas
with handle connector



S26-202-01
1,7 mm Ø
5 cm

S26-203-01
1,7 mm Ø
15 cm



S26-204-02
2 mm Ø
5 cm

S26-205-02
2 mm Ø
15 cm



S26-206-03
3 mm Ø
5 cm

S26-207-03
3 mm Ø
15 cm



S26-208-01
1,7 mm Ø
5 cm

S26-209-01
1,7 mm Ø
15 cm

S26-210-02
2 mm Ø
5 cm

S26-211-02
2 mm Ø
15 cm

S26-212-03
3 mm Ø
5 cm

S26-213-03
3 mm Ø
15 cm

Cannulas
with Luer-Lock connection

Measurements = working length

Handles with connection thread for cannulas and hose connection 0 13 mm, for operation with a suction pump.

Recommended handles for the tumescence local anaesthesia technique
Hole for interrupting the suction process



S26-214-00
Fix Handle



S26-215-00
with rotatable tube connection
Handle



S26-216-00
Fix Handle



S26-217-00
with rotatable tube connection
Handle

Liposuction Cannulas
Cánulas de liposucción
Liposuction-Kanülen
Canules pour liposucccion

Liposuction cannulas, stainless steel with thumb plate and hose connection Ø 13 mm, for use with a suction pump.



S26-221-03
 3 mm Ø
 15 cm
 one central hole

S26-222-04
 4 mm Ø
 25 cm
 one central hole

S26-223-06
 6 mm Ø
 15 cm
 one central hole

S26-224-06
 6 mm Ø
 25 cm
 one central hole

S26-225-09
 9 mm Ø
 30 cm
 one central hole

**Gasparotti
 Cannulas**

Measurements = working length



S26-226-04
 4 mm Ø
 15 cm
 two central holes

S26-227-04
 4 mm Ø
 25 cm
 two central holes

S26-228-06
 6 mm Ø
 30 cm
 two central holes

**Gasparotti
 Cannulas**

Measurements = working length



S26-229-03
 3 mm Ø
 15 cm
 one central hole

S26-230-06
 6 mm Ø
 15 cm
 one central hole

S26-231-06
 6 mm Ø
 25 cm
 one central hole

**Gasparotti
 Cannulas**

Measurements = working length

Liposuction Cannulas Set

S26-235-00

Liposuction Cannula Set for Abdomen and Saddlebag With Threaded Fitting
 Liposuction Cannula with threaded connector, 22-Holes, diameter 3.0 mm, working length 20cm

Liposuction Cannula with threaded connector, 22-Holes, diameter 4.0 mm, working length 20cm

Liposuction Cannula with threaded connector, 30-Holes, diameter 3.0 mm, working length 30cm

Liposuction Cannula with threaded connector, 30-Holes, diameter 4.0 mm, working length 30cm

Liposuction Cannula Handle

Liposuction Cannula Handle with Interrupting Hole

S26-240-00

Liposuction Cannula Set for Face with Luer Lock Fitting

Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 5 cm

Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 15 cm

Liposuction Cannula, 18-Holes, diameter 2.0 mm, working length 15 cm


S26-245-00

Liposuction Cannula Set for Face with Threaded Fitting

Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 5 cm

Liposuction Cannula, with one central hole, diameter 1.7 mm, working length 15 cm

Liposuction Cannula, 18-Holes, diameter 2.0 mm, working length 15 cm

Liposuction Cannula Handle

Liposuction Cannula Handle with Interrupting Hole

S26-250-00

Liposuction Cannula Set for Neck

Liposuction Cannula with Luer-Lock, with one central hole, diameter 2.0 mm, working length 5cm

Liposuction Cannula with Luer-Lock, with one central hole, diameter 3.0 mm, working length 5cm

Liposuction Cannula with Luer-Lock, with one central hole, diameter 3.0 mm, working length 15cm

Liposuction Cannula with threaded connector, 18-Holes, diameter 3.0 mm, working length 15 cm

Liposuction Cannula with threaded connector, Angled, 22-Holes, diameter 3.0 mm, working length 20 cm

Liposuction Cannula with threaded connector, 8-Holes, diameter 2.0 mm, working length 15 cm

Liposuction Cannula Handle

Liposuction Cannula Handle with Interrupting Hole

Liposuction Cannulas Set



S26-255-00

Liposuction Cannula Set for Thie and Lower Leg with Threaded Fitting

Liposuction Cannula with threaded connector, byonet shape, 22-Holes, diameter 3.0 mm, working length 20 cm

Liposuction Cannula with threaded connector, byonet shape, 22-Holes, diameter 4.0 mm, working length 20 cm

Liposuction Cannula with threaded connector, byonet shape, 30-Holes, diameter 3.0 mm, working length 30 cm

Liposuction Cannula with threaded connector, byonet shape, 30-Holes, diameter 4.0 mm, working length 30 cm

Liposuction Cannula Handle

Liposuction Cannula Handle with Interrupting Hole

S26-260-00

Micro Fat Injection Cannulas Set

Open Port Injector Cannula, 1.2mm, 10cm, Luer Lock

Multi Holes Harvester Cannula Model 2, 2.5mm, 20cm, Luer Lock

Tumescent Infiltrator Cannula, 1.5mm, 12cm, Luer Lock

Micro Injector Cannula, 0.8mm, 4cm, Luer Lock

Micro Injector Cannula, 1.0mm, 5cm, Luer Lock

Holding Device

Luer Lock Syringe Cap

Safe Cannula Set

S26-265-00

1 TWO STANDARD HOLE LIPOSUCTION CANNULA, 3.5mm, 25cm, Super Handle

2 Mercedes Liposuction Cannula, 5.0mm, 30cm, Super Handle

3 Mercedes Liposuction Cannula, 4.5mm, 30cm, Super Handle

4 Mercedes Liposuction Cannula, 4.0mm, 25cm, Super Handle

5 Single Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle

6 Double Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle

7 Triple Basket Liposuction Cannulas, 4.0mm, 30cm, Super Handle

8 Spiral Three Port Harvester Cannula, 3.0mm, 20cm, Super Handle

9 Spiral Three Port Harvester Cannula, 4.0mm, 25cm, Super Handle

10 Three Port Harvester Cannula, Model 4, 4.5mm, 30cm, Super Handle

11 Three Ports Pyramid Cannula, 3.5mm, 25cm, Super Handle

12 Basket Cannula, 5.0mm, 32cm, Super Handle




Fettinjektions- und Fettabsaugspritzen

Jeringas para inyección de grasa y succión de grasa

Seringues d'injection de graisse et d'aspiration de graisse

Siringhe per iniezione e aspirazione di grasso

Fat Injection Syringes




S26-301-50
Fat Injection Syringe
50ml



S26-302-20
Universal Fat Injection syringe
20ml



S26-303-10
Universal Fat Injection syringe
10ml




S26-304-01
Fat Injection Syringes
0.1 ml per click
1 ml



S26-305-03
Fat Injection Syringe
3 ml special
0.1 ml per click

Fat Suction Syringes




S26-311-50
Fat Suction Syringe
50 ml




S26-312-50
Fat Suction Syringe
50 ml



S26-313-20
Fat Suction Syringe
20 ml



S26-314-10
Fat Suction Syringe
10 ml



S26-315-05
Fat Suction Syringe
5 ml



S26-316-03
Fat Suction Syringe
3 ml



LEADER SURGICAL

Manufacturers & Exporters of Surgical

42 A-1 Gulberg 3 Lahore Pakistan

Cell: + 92-322-3000333

Email: info@leadersurgical.org

Website: www.leadersurgical.org

