



**DIAGNOSTIC  
INSTRUMENTS**



S1-101-00



S1-102-00



S1-103-00



S1-104-00



S1-105-00



S1-106-00



S1-107-00

#### Otoscope Fiber Optic

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting.  
Fixed lens 6X in plastic mount may be removed for introduction of instruments.  
Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Ionized  
Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v  
Available in multi colours.

**Colour Choice**

Red	
Green	
Black	
Blue	
Yellow	
Orange	



S1-108-00  
Mini Ophthalmoscope



S1-116-00  
Mini Ophthalmoscope





**S1-117-00**

Varioflex combination hammer according to Fassbender, with never needle and brush. Solid handle with centimeter scale marking (approx. 155g) length 20 cm.



**S1-118-00**

Rossier - Reflex hammer, with elastic handle, length: 36cm (approx. 95g)



**S1-119-00**

Rossier - Reflex hammer, with elastic handle, length: 24cm (approx. 90g)



**S1-120-00**

Babinski - Reflex and percussion hammer, interchangeable with screw fitting, with metal grip (approx. 95g)



**S1-121-00**

Babinski - Reflex and percussion hammer, interchangeable with screw fitting, with metal grip (approx. 90g)



**A**

**S1-122-00**



**B**

**S1-123-00**



**C**

**S1-124-00**



**D**

**S1-125-00**

Colorflex, reflex and percussion hammer, Small version, length: 17cm  
Color: red, blue, yellow, black (approx. 100g)



**DUPLEX**  
**Stethoscope**  
**S1-137-00**  
super lightweight



**DUPLEX-BABY**  
**Stethoscope**  
**S1-138-00**  
super lightweight



**ANAESTOPHON**  
**Stethoscope**  
**S1-139-00**  
super lightweight



**Clar Mead Mirror**  
**komplett**  
**6 Volt**  
**S1-151-00**



**Ziegler**  
**S1-152-00**  
**Ø 90mm**  
with head band without  
protection cap



**Ziegler**  
**S1-153-00**  
**Ø 90mm**  
with head band with  
protection cap



**S1-156-00**  
**Ø 50mm**



**PINARD WOODEN**  
**Stethoscope**  
**S1-157-00**  
18 cm



**PINARD ALUMINIUM**  
**Stethoscope**  
**S1-158-00**  
15 cm



**WEBER**  
**Aesthesiometer**  
**S1-161-00**  
9 cm



**MAGILL Catheter Intubating Forcep GOLD**

**S1-162-20** 20cm (6.5cm shaft), Infant, hole (dia 2.00 mm)

**S1-163-23** 23cm (7.5cm shaft), Infant, hole (dia 2.30 mm)

**S1-164-26** 26cm (10cm shaft), Infant, hole (dia 2.60 mm)



**MAGILL**  
Catheter Introducing  
Forcep  
**S1-191-17** 17cm



**MAGILL**  
Catheter Introducing  
Forcep  
**S1-192-20** 20cm




**MAGILL**  
Catheter Introducing  
Forcep  
**S1-193-25** 25cm



**NEURO-DUOFLEX**  
Hammer  
S1-201-22 22cm



**DEJERINE**  
Hammer  
S1-202-20  
20cm



**DEJERINE**  
Hammer  
S1-203-20  
20cm



**TROEMNER**  
Hammer  
S1-204-14 145 g  
S1-204-17 175 g  
S1-204-23 235 g  
24.5cm  
S1-204-00 spare rubber



**TROEMNER-SIMPLEX**  
Hammer  
S1-205-21  
21cm



**BUCK Hammer**  
S1-206-18 18.5cm  
S1-206-20 20cm



**BUCK HAMMER**  
S1-207-18 18cm  
S1-207-20 20cm



**TAYLOR Hammer**  
S1-208-16 16cm  
S1-208-18 18cm  
S1-208-20 20cm



**TAYLOR Hammer**  
S1-209-18  
18cm



Schlaghämmer, Windrad

Martillos de percusión, molinete

Marteaux à percussion, moulinet

Martelli a percussione, girandola



**Babinsky Hammer**  
S1-216-22 22cm



**Babinsky Hammer**  
S1-217-22 22cm




**Rabiner Hammer**  
S1-218-25 25.5cm  
with needle & brush



**National Hospital type  
Percussion Hammers**  
Length diameter  
S1-219-50 350mm, 50mm  
S1-219-60 350mm, 60mm  
S1-219-70 350mm, 70mm



**Berliner Hammer**  
S1-226-20 20cm




**Traube Hammer**  
S1-227-16 16cm



**Wartenberg Pinwheel**  
S1-231-19 19cm

## Conventional Laryngoscope Set

### MC-INTOSH LARYNGOSCOPE



#### BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

#### Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel Matt Type, Standard Illumination 2.5 V.

#### Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

#### Mc-INTOSH LARYNGOSCOPE SET STANDARD ILLUMINATION

##### S1-300-00

##### PEDIATRIC SET

Set of 3 Blades

Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

##### S1-301-00

##### STANDARD SET

Set of 3 Blades

Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

##### S1-302-00

##### ADULT SET

Set of 4 Blades

Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

##### S1-303-00

##### ADULT SET

Set of 5 Blades

Fig. 0,1,2,3,4  
Supplied in plastic case with 2 spare blade

Fig. 0

S1-304-00

Fig. 1

S1-304-01

Fig. 2

S1-304-02

Fig. 3

S1-304-03

Fig. 4

S1-304-04



S1-304-20  
Spare Bulb




S1-304-10  
Laryngoscope Battery  
Handle for Adult  
C type 2 Batteries



Mc-Intosh Laryngoscope Set  
Standard Illumination  
S1-302-00

## Conventional Laryngoscope Set

### MILLER LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

#### Mc-Intosh Laryngoscope

Plastic case, Standard illumination.

Manufactured from AISI 304 Stainless Steel.

#### Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

### MILLER LARYNGOSCOPE SET OF CONVENTIONAL STANDARD ILLUMINATION

#### S1-305-00 PEDIATRIC SET

Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

#### S1-306-00 STANDARD SET

Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

#### S1-307-00 ADULT SET

Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

#### S1-308-00 ADULT SET

Set of 5 Blades  
Fig. 0, 1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 0  
S1-309-00



Fig. 1  
S1-309-01



Fig. 2  
S1-309-02



Fig. 3  
S1-309-03



Fig. 4  
S1-309-04



S1-309-20  
Spare Bulb




S1-309-10  
Laryngoscope Battery  
Handle for Adult C type  
2 batteries



Laryngoscope Set of conventional Miller  
Standard Illumination  
S1-307-00



**MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE**

**BLADE**

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

**Mc-Intosh Laryngoscope Fiber Optic**

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades.

**Blades Availability:**

- \* Matt Finish
- \* Mirror Finish
- \* Sand Blast Finish
- \* Pls, confirm in your order list

**Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION**
**S1-320-00**
**PEDIATRIC SET**

Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade

**S1-321-00**
**STANDARD SET**

Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade

**S1-322-00**
**ADULT SET**

Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade

**S1-323-00**
**ADULT SET**

Set of 5 Blades  
Fig. 0,1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 0  
S1-324-00



Fig. 1  
S1-324-01



Fig. 2  
S1-324-02



Fig. 3  
S1-324-03



Fig. 4  
S1-324-04



S1-324-20  
Laryngoscope  
Bulb Small



S1-324-22  
Laryngoscope  
Bulb Large




S1-324-10  
Laryngoscope Battery  
Handle for Adult  
C type 2 batteries



S1-322-00  
Mc-Intosh Fiber Optic Laryngoscope Set  
Halogen Illumination



**MILLER FIBER OPTIC LARYNGOSCOPE**

**BLADE**

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

**MILLER FIBER OPTIC LARYNGOSCOPE**

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades.

**Blades Availability:**

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

**MILLER FIBER OPTIC LARYNGOSCOPE HALOGEN ILLUMINATION**

**S1-325-00**  
**PEDIATRIC SET**  
Set of 3 Blades  
Fig. 0,1,2  
Supplied in plastic case with 2 spare blade



Fig. 0  
**S1-329-00**



Fig. 1  
**S1-329-01**

**S1-326-00**  
**STANDARD SET**  
Set of 3 Blades  
Fig. 1,2,3  
Supplied in plastic case with 2 spare blade



Fig. 2  
**S1-329-02**

**S1-327-00**  
**ADULT SET**  
Set of 4 Blades  
Fig. 1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 3  
**S1-329-03**

**S1-328-00**  
**ADULT SET**  
Set of 5 Blades  
Fig. 0,1,2,3,4  
Supplied in plastic case with 2 spare blade



Fig. 4  
**S1-329-04**



**S1-329-20**  
Laryngoscope Bulb Small



**S1-329-22**  
Laryngoscope Bulb Large



**S1-329-10**  
Laryngoscope Battery  
Handle for Adult  
C type 2 Batteries



**S1-328-00**  
**MILLER FIBER OPTIC LARYNGOSCOPE**

## Laryngoscope Flexible (Mc-Coy)

### CONVENTIONAL SET

**S1-351-00 Flexible Laryngoscope Conventional Set of 3 blades**  
 Fig. 2,3,4

**S1-352-00 PEDIATRIC SET**  
 Handle With Blade Fig. 2  
 Supplied in plastic case

**S1-353-00 STANDARD SET**  
 Handle With Blade Fig. 3  
 Supplied in plastic case

**S1-354-00 ADULT SET**  
 Handle With Blade Fig. 4  
 Supplied in plastic case  
 Halogen Illumination  
 Manufactured From  
 Stainless Steel High  
 grade 304 material.



**S1-355-10 'C' type**  
 battery handle  
 Blades Size:  
**S1-355-02** Fig. 2  
**S1-355-03** Fig. 3  
**S1-355-04** Fig. 4

LEADER  
 SURGICAL

**S1-361-00 Flexible Laryngoscope Conventional Set of 3 blades**  
 Fig. 2,3,4

**S1-362-00 PEDIATRIC SET**  
 Handle With Blade Fig. 2  
 Supplied in plastic case

**S1-363-00 STANDARD SET**  
 Handle With Blade Fig. 3  
 Supplied in plastic case

**S1-364-00 ADULT SET**  
 Handle With Blade Fig. 4  
 Supplied in plastic case

Halogen Illumination  
 Manufactured From  
 Stainless Steel High  
 grade 304 material  
 packed in green plastic  
 box.



**S1-364-00**



**S1-365-10 'AA' type**  
 battery handle  
 Blades Size:  
**S1-365-02** Fig. 2  
**S1-365-03** Fig. 3  
**S1-365-04** Fig. 4



Two laryngeal  
mirrors: n. 4 and 5



Handle

Bent-arm illuminator:  
pin contact fitting.

Specula



With 2 spare Bulbs

**Professional Otoscope & Ophthalmology Set**  
**ENT Medical Diagnostic Set**  
**S1-400-00**

2.5 V Standard illumination for use with  
2 batteries C type.

Including:

Otoscope Head: Luer Lock fitting.

Chrome plated.

Fix lens 4 x in plastic mount can be  
removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luer Lock fitting.

Lens selection wheel permits choice of corrective  
lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle aperture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in Quality Plastic Case.



Otoscope Head



Nasal speculum: Chrome plated.



Ophthalmoscope head



ENT Specula





**Universal Otoscope & Ophthalmology Set**  
**ENT Medical Diagnostic Set**  
**S1-410-00**

2.5 V Standard illumination for use with  
 2 batteries type C.

Including:

Otoscope Head: Luer Lock fitting.  
 Chrome plated.

Fix lens 4 x in plastic mount can be  
 removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: pin contact fitting.

Lens selection wheel permits choice of corrective  
 lenses from 0 to - 20 and 0 to +20 diopters.

Single large circle operture.

Bent-arm illuminator: Luer Lock fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in Quality Plastic Case.



Two laryngeal  
 mirrors: n. 4 and 5



Bent-arm illuminator:  
 pin contact fitting.



Nasal speculum: Chrome plated.



ENT Specula



Otoscope Head



Ophthalmoscope Head



Handle




**DIAGNOSTIC PENTOSCOPE SET  
S1-420-00**

2.5 V Standard illumination for use with 2 batteries type AA.

Consisting:

Universal Handle thread type fit for all Pentoscope Heads like.

- Otoscope Head.
- Ophthalmoscope Head.
- Throat Branch Head.
- Mouth Mirrors Head.

It have two exploring mirrors No. 4 & 5 mm.

Its spare 06 disposable Plastic ear Canulas are with sizes 2 x 2.50 mm, 2 x 3.50 mm & 2 x 4.5 mm.

On special order it can be as 2 x 1.5 mm for new born babies 2 x 2.5 mm for pediatric and 2 x 4 mm for adults, so it means total 06 canula also available in a plastic tube packing.

Sturdy adjustment mount bracket

2 spare 2.5 V bulb.

Pentoscope set battery handle is operating with two AA size batteries.

This Pentoscope set Packed in our new beautiful Blue Colour plastic case in same blue inner packing material.

Spare parts also available of any item.



With 2 spare Bulbs



Mouth Mirrors Head



Standard illumination

LEADER  
SURGICAL



Otoscope Head



Sturdy Adjustment Mount Bracket



Plastic ear Canulas

*Diagnostik-Bestestecke*
*Juegos para diagnóstico*
*Trousses à diagnose*
*Astucci diagnostici*


**Standard Otoscope  
S1-425-00**

This otoscope is 2.5 V standard illumination, handle is powder coated as well otoscope head too. Handle and head is Pin Screw lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Nice Wooden Box Packing



**Standard Otoscope  
S1-435-00**

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Standard Otoscope  
S1-430-00**

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle is Pin Screw lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Mini Otoscope Fiber Optic  
S1-440-00**

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Ionized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colour.



**Mini Otoscope Fiber Optic  
S1-445-00**

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Ionized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colour.







**Standard Otoscope + Ophthalmoscope Set  
S1-450-00**

This Set is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Box Packing.



Otoscope Head



ENT Specula



Ophthalmoscope head



ENT Specula



'C' type Battery handle



**Standard Ophthalmoscope Set  
S1-455-00**

A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Packed in Black Plastic Case.



Ophthalmoscope head



'C' type Battery Handle


**Standard Otoscope Set  
S1-460-00**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing


**Standard Otoscope Set  
S1-465-00**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Reusable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



**Otoscope + Ophthalmoscope Set  
S1-470-00**

This otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Raisable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing.




**Mini Ophthalmoscope  
S1-475-00**

A Full Functional Ophthalmoscope.  
compact durable, light weight.  
Standard bulb and precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
Made of plastic Material.

**Colour Choice**

**Mini Ophthalmoscope  
S1-480-00**

A Full Functional Ophthalmoscope.  
compact durable, light weight.  
Standard bulb and precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
made of plastic Material.


**Standard Ophthalmoscope  
S1-485-00**

A Full Functional Ophthalmoscope.  
compact durable, Standard bulb and  
precision system  
4 different apertures selection  
range of lenses -20 to +20  
bright white light and brilliant image  
Handle made of Stainless Steel Material.

*Messinstrumente*
*Instrumentos de medidas*
*Instruments pour mesurer*
*Strumenti di misura*

**GRADUATED TAPE MEASURER**
**S1-501-02**


 2 m long (1.5 m) graduated  
 in inch & cm markings

**S1-503-15** 15cm  
**S1-503-30** 30cm  
 cm & mm  
 markings

**JAMESON CALIPER**  
**S1-505-08**  
 8cm

**TOWNLEY CALIPER**  
**S1-507-10**  
 10cm

LEADER  
SURGICAL


**VERNIER CALIPER**  
**S1-509-15**  
 15cm

**VERNIER CALIPER**  
**S1-511-15**  
 15cm

**MÖELTGEN Flexometer**  
**S1-517-24**  
 24cm

**CARROLL POCKET GONIOMETER**  
**S1-519-14**  
 14.5cm

**BERTILLON MEASURING DEVICE**  
**S1-521-32**  
 32cm



**CASTROVIEJO-EPKER**  
**Caliper**

**S1-525-09** 9 cm Str.  
 Measuring range, 40 mm  
**S1-526-09** 9 cm Cvd.  
 Measuring range, 40 mm



**CASTROVIEJO-EPKER**  
**Caliper**

**S1-528-18** 18 cm  
 Measuring range  
 40 mm, straight



**CASTROVIEJO-EPKER**  
**Caliper**

**S1-530-18** 18 cm Angled  
 Measuring range  
 40 mm



**CASTROVIEJO**  
**S1-532-18** 18cm  
 Bone Caliper  
 Straight



**CASTROVIEJO**  
**S1-534-18** 18cm  
 Bone Caliper  
 Curved



**CASTROVIEJO**  
**S1-539-08**  
 8 cm  
 Eye View Caliper  
 Measuring range  
 20 mm



**Tonometer**  
**S1-545-00**  
 Set





# LEADER SURGICAL

**Manufacturers & Exporters of Surgical**

42 A-1 Gulberg 3 Lahore Pakistan

Cell: + 92-322-3000333

Email: [info@leadersurgical.org](mailto:info@leadersurgical.org)

Website: [www.leadersurgical.org](http://www.leadersurgical.org)

