





**DISSECTING KITS
LABORATORY
INSTRUMENTS**


S28-101-00

Dissecting Set of 9 Pcs.

- Dissecting Forceps 11.5cm Sharp Point.
- Teasing needle Sharp STR.
- Teasing Needle Sharp Angular.
- B.P Handle No.4.
- Blade No. 22.
- Thumb Forcep T.O.E 13cm.
- Dissecting Scissor STR 11.5cm.
- Dissecting Scissor STR 13cm.
- Zipper Case Inside Red Velvet.


S28-102-00

Dissecting Set of 10 Pcs.

- Grooved Director with Tongue Tie 14cm.
- Probe Double 14cm.
- Mayo Hegar Needle Holder 14cm.
- B.P Handle No.3
- B.P Handle No.4
- Dressing Scissor 14cm S/B STR.
- Dressing Scissor 14cm S/B CVD.
- Dressing Forcep 1x2 teeth 14cm.
- Dressing Forcep Serrated 14cm.
- Zipper Case Inside Red Velvet.


S28-103-00

Dissecting Set of 10 Pcs.


- Probe Double 13cm.
- 2PCS. Scalpel S.S.
- Grooved Director With Tongue Tie 14cm
- Mayo Hegar Needle Holder 14cm.
- Dissecting Forceps 14cm 1x2 Teeth.
- Dressing Forceps 14cm Serrated.
- Dressing Scissor 14cm S/B STR.
- Dressing Scissor 14cm S/B CVD.
- Zipper Case Inside Red Velvet.


S28-104-00


Dissecting Set of 10 Pcs.

- Blow Pipe.
- 3PCS Scalpel.
- Dissecting Forceps 14cm.
- Grooved Director With Tongue Tie 14cm
- Triple Chain Hook.
- Dressing Scissor Probe Pointed 14cm.
- Probe Double 14 cm.
- Zipper Case Inside Red Velvet.


Dissecting Set
S28-105-00 Containing 9 Instruments


S28-106-00
Basic Dissecting Set

A Collection of Instruments For Introductory Dissection vinyl pouch
 Scalpel 6" Scissors 4 1/2" Forceps student grade 4 1/2"
 Straight needle in Wooden Handle Bent needle in wooden handle
 Medicine Droppers.


S28-107-00
Dissecting Set

Beginners For Elementary Junior High Or High School Life Science
 Vinyl pouch Scissors Forceps student /Scalpel with Replaceable
 blade Straight needle in wooden handle Bent needle in wooden
 handle Medicine dropper


S28-108-00
General Biology Dissecting Set

Designed for work in general biology good Quality instruments
 scissor student Forceps fine point 4 1/2" S/S
 Scalpel 1 1/2 Blade (with G pin S/S)
 Straight Needle with wooden handle Medicine Dropper
 single fold leatherette case.


S28-109-00
Fine Dissecting Set

Designed for Freshman Courses in Biology
 Scissors Fine point 4 1/2" Forceps fine point with g pin 4 1/2"
 Scalpel replaceable Blade Teasing Needle Str
 Teasing Needle Curved Medicine Dropper
 Single Fold Leatherette Case.


S28-110-00
Intermediate Dissecting Set

Contains high quality instruments recommended for advance
 biology courses . Scissor surgical 4 1/2" Forceps fine point 4 1/2 "
 Surgeon scalpels (with # 22 blade) Teasing needle plastic handle Str
 Teasing needle curved Plastic handle Mall probe and seeker
 Medicine dropper Single fold leatherette case


S28-111-00
Tissue Culture Dissecting Set

Features Stainless Steel Instrument in leatherette case
 Scissors M/Dissecting Iris 4 1/2 " Forceps watchmaker 4 1/4"Str
 Forceps watch make 4 1/4" curved Scalpel micro dissecting 5"
 18x3 mm blade Teasing Needle Str S/S Teasing Needle curved S/S
 Single Fold leatherette case


S28-112-00
Anatomy Dissecting set

A complete set Detailed Dissection for more advanced work in zoology and comparative anatomy .
 Surgical scissors 5 1/2". Fine point scissors 4 1/2"
 Forceps fine point curved 4 1/2" Thumb forceps superior 5"
 Scalpel 1 1/2" Blade .
 Cartilage knife/Flexible Probe. Triple chain Hook/Mall probe.
 Needle wooden handle Str. Needle wooden handle Str
 Medicine Dropper. Double fold leatherette case


S28-113-00
Graduate Dissecting Set

Design for the biology graduate students
 Surgical scissors 5 1/2". Fine point scissors 4 1/2".
 Forceps fine point Str 4 1/2" Forceps fine point curved 4 1/2".
 Tripple chain hook. cartilage knife scalpel stain finish 1 1/2" blade
 forceps cartilage 5" 2 needle holder with straight needle .
 Flexible probes. Probes with hook. Probe & seeker
 Dropping Pipe Double Fold Leatherette case


S28-114-00
Micro Dissecting Set

Feature stainless steel Instruments in Leatherette case.
 Scissors M/Dissecting Iris 4 1/2" Scissors sharp point Iris 4 1/2"
 Forceps Micro Dissecting 4 1/4" Straight Forceps micro dissecting 4"
 Curved Scalpel micro dissecting 18 x 3 mm blade
 scalpel micro dissecting 30 x 3 mm blade Dissecting needle Str .
 Dissecting needle Curved. Dissecting needle Full-Spear point
 Dissecting needle Half-Spear point In double folding case.


S28-115-00
Large Animal set

This set provide all the Necessary instruments Dissection of large Animal such as Cat, Shark and Rabbits Scissor surgical 5 1/2" D/Scissors superior 4 1/2" Scissors fine point Str 4 1/2" Forceps large Tissue 1x2 teeth 6" Forceps utility curved point 6" .
 Forceps fine point str 4 1/2" Scalpel screw Lock with # 22 blade .
 Knife cartilage 6 1/2". Probe sharp point angular. Section lifter S/S
 Teasing needle Str Wooden handle. Teasing Needle Curved Wooden handle
 Medicine dropper.


S28-116-00
New Instructor's Dissecting set

Contain quality stainless steel Instrument Enclosed in an attractive zippered case Scissor surgical 4 1/2" scissors fine point 4 1/2" Forceps fine point Str. Forcep fine point Curved. Mosquito forceps Curved 5" box lock. One piece scalpel 1 1/2" blade Surgeon handle # 4. Mall probe and seeker 2 plastic handle str needle also 2 blade # 22 for each set.


S28-117-00
Dissecting Set - Advanced

Stainless steel, consisting of 12 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Artery Forceps - 1 pc. Tissue Forceps - 1 pc. Seeker - 1 pc. Scalpel Handle w/ Blade - 1 pc. Bone Cutter - 1 pc. Dental Probe & Seeker - 1 pc


S28-118-00
Dissecting Kit

Kit includes:

1 4 1/2" Dissecting forceps with curved and serrated tips
 1 Curved dissecting needle 1 Glass pipette with rubber bulb
 1 4 1/2" Dissecting forceps straight and serrated tips
 1 Dissecting scalpel, 6" long, 1 1/2" cutting edge.
 1 4 1/2" Dissecting scissors, sharp, straight points
 1 Plastic ruler 6", 1/16" div., millimeter and centimeter


S28-119-00
Dissecting Set - Pro

Stainless steel, consisting of 20 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc. B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Ruler 6" - 1 pc. Magnifying Glass - 1 pc. Glass Slide - 1 pc. Chain and Hook - 1 pc. Cover Slips - 1 pc. Blow Pipe - 1 pc. Needle Wooden Handle - 1 pc. Dropper Glass - 1 pc. Brush - 1 pc.


S28-120-00
Dissecting Set - Popular

Stainless steel, consisting of 14 instruments for accurate dissection. Complete in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needles - 4 pcs. Scalpel - 1 pc. B.P. Handle w/ Blade - 1 pc. Razor - 1 pc. Magnifying Glass - 1 pc. Chain and Hook - 1 pc. Blow Pipe - 1 pc.


S28-121-00
Dissecting Set - Senior

Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

Scissors - 2 pcs. Forceps - 2 pcs. Needle - 2 pcs. Scalpel - 1 pc. Scalpel Handle w/ Blade - 1 pc. Ruler 6" - 1 pc.

Diagnostik-bestecke

Diagnostic Setts

Estuches de diagnostico

Combinazioni per diagnostica


S28-122-00
Dissecting Set - Student

Stainless steel, consisting of 9 instruments for accurate dissection. Supplied in zipper leather case.

Kit consists of:

 Scissors - 2 pcs. Forceps - 2 pcs.
 Needle Wooden Handle- 2 pcs.
 Scalpel - 1 pc. Scalpel Handle w/ Blade - 1 pc.
 Dropper Glass - 1 pc.

S28-123-00
Dissecting Set - Beginner

Stainless steel, consisting of 7 instruments for accurate dissection. Complete in zipper leather case.

Kit consists of: Scissors - 2 pcs.

Forceps - 2 pcs. Seeker - 1 pc.

B.P. Handle w/ Blade - 1 pc.

Needle Wooden Handle - 1pc


S28-124-00
MICROSCOPIC DISSECTING SET

S28-125-00
Dissecting-kits
11 Pcs

S28-126-00
Biology Dissecting Tools Set

S28-127-00
Anatomical Dissecting Set

 13 parts filled, in leatherette zipper pocket
 Content

1 straight dissecting needle with plastic handle, 140 mm

1 bent dissecting needle with plastic handle, 140 mm

1 lancet with plastic handle, 140 mm

1 bandage scissors, pointed/blunt, 130 mm

1 straight microscope scissors, 115 mm

1 bent microscope scissors, 105 mm

1 straight pointed tweezers, 105 mm

1 straight blunt tweezers, 115 mm

1 pointed scalpel, 130 mm, cutting edge 17 mm

1 convex scalpel, 130 mm, cutting edge 17 mm

1 pointed scalpel, 155 mm, cutting edge 40 mm

1 section lifter, 160 mm, width 10 mm

1 directional probe, 140 mm


S28-128-00
DE-SET1


S28-129-00
Dissecting Set - Economy

Stainless steel, consisting of 7 instruments for accurate dissection. Supplied in leather storage case.
 Kit consists of: Scissors - 1 pc. Forceps - 1 pc.
 Needles - 2 pcs. B.P. Handle w/ Blade - 1 pc.
 Dropper Glass - 1pc. Ruler 6" - 1 pc.


S28-130-00
Basic Dissecting Kit

S28-131-00
Basic Dissecting Tools Set

S28-132-00
Deluxe Instructor's Dissecting Set

S28-133-00
Podiatry Kit

S28-134-00
Anatomy Dissecting Set

S28-135-00
Micro Dissecting Set

Features stainless-steel instruments in Double Leatherette Case



S28-136-00
Standard Dissecting Set II



S28-137-00

Thomas Basic Dissecting Kit

Designed for students or limited use 6 instruments
 Chrome and stainless steel construction Kit includes:

6" chrome scalpel

4 1/2" Stainless steel dissecting scissors

4 1/2" Chrome-plated forceps

Straight dissecting needle Bent dissecting needle

Pipette 6" Ruler Single fold vinyl case with lining



S28-138-00
General Biology Dissecting Set



S28-139-00
Dissecting Set, Biology Training, Molded Plastic Case



S28-140-00
Advanced Anatomy Dissecting Set



S28-141-00
**Advanced Anatomy Dissecting Set
 with Leather Case**

Diagnostik-bestecke

Estuches de diagnostico

Diagnostic Setts

Combinazioni per diagnostica



S28-142-00
Teacher's Dissecting Set



S28-143-00 Comparative Anatomy Dissecting Sets
Set is sold with either Double Leatherette Case or Molded Plastic Case.



S28-144-00
Dissecting Set



S28-145-00
MICROSCOPIC DISSECTING SET
5 parts or 7 parts filled, in wooden case Content
1 straight scissors, 100 mm
1 bent dissecting needle, 150 mm
1 straight dissecting needle, 150 mm
1 tweezers, straight point, with guide pin, 115 mm
1 scalpel 125 mm



S28-146-00
CHEMISTS' TOOL SET
4 parts filled, in a case with press stud Content
1 powder spatula, 150/9 mm
1 micro-powder spatula 130/4 mm
1 pointed/blunt scissors, 130 mm
1 straight pointed tweezers, 130 mm



S28-147-00
Dissecting Kit



S28-151-00
Dissecting Set



S28-152-00
Dissecting Set



S28-153-00
Dissecting Set



S28-154-00
Dissecting Set



S28-155-00
Dissecting Set



S28-156-00
Dissecting Set



S28-157-00
Dissecting Set



S28-158-00
Dissecting Set



S28-159-00
Dissecting Set



S28-160-00
Dissecting Set



S28-161-00
Dissecting Set



S28-162-00
MICROSCOPIC DISSECTING SET



S28-163-00
Deluxe Dissecting Tools Set



S28-164-00
Dissecting Set



S28-165-00
DISSECTING SET



S28-166-00
DISSECTING SET




S28-167-00
DISSECTING SET



S28-168-00
DISSECTING SET







S28-179-00
Dissection Kit w/ Screw Lock Blade




S28-180-00
Precision Dissection kit
Precision Dissection kit - plastic case



S28-181-00
Dissection Kit - Basic



S28-182-00
Budget Dissection Kit

	 S28-183-00 Biology Dissection Kit
	 S28-184-00 Anatomy Dissection Kit & Stethoscope Package Med-Student first year gross anatomy dissection kit
	 S28-185-00 Zippy Precision Kit
	 S28-186-00 Introductory Anatomy Dissection Kit

	 S28-187-00 Fine Zippy Dissection Kit
	 S28-188-00 Precision Plus Kit
	 S28-189-00 Zippy Dissection Kit - 65ZP
	 S28-190-00 Dissection Kit - Intermediate II - Kit-1PC



S28-191-00
Precision Dissection Kit Package with
Safety Eyewear and Gloves



S28-192-00
Comprehensive Dissection Kit w/
Goggle - Kit-5



S28-193-00
Entomology Dissection Kit
for insect studies.



S28-194-00
Dissection Kit - Intermediate III - Kit-2

	 S28-195-00 Student Dissection Kit (Plastic Case)
	 S28-196-00 Zippy Dissection Kit w/ Screw Lock Scalpel
	 S28-197-00 Veterinary Student Dissection Kit First Year Veterinary Student Anatomy Training Kit
	 S28-198-00 Podiatry Students Training Kit Podiatrist Student Dissecting Kit


S28-199-00
Zippy Dissection Kit Package With Nitril Gloves and Safety Eyewear

Ideal dissection package includes Zippy dissection set, eyewear, and pack of nitrile gloves


S28-200-00
Comprehensive Economy Dissection Kit

Comprehensive Dissecting kit includes 13 instruments and 2 packs of blades to set any student up for a science lab. All tools come neatly packed away in a snap button case.


S28-201-00
Student Anatomy Dissection Kit - Economy

S28-202-00
Invertebrate Dissection Kit

Includes 8 tools and 2 blades in a sleek vinyl zippered case. All instruments are stainless steel unless otherwise noted.

	 S28-203-00 Budget Entomology Kit Budget Entomology dissection kit includes the basic tools needed for a classroom dissection.
	 S28-204-00 General Biology Dissection. Kit Basic General Biology Dissecting Kit
	 S28-205-00 Med Student Anatomy Dissection Kit Anatomy dissection kit for first year medical students!
	 S28-206-00 Dissection Kit Standard



S28-207-00
Med Student Anatomy Dissection Kit
Anatomy Dissection Kit for
first year medical students



S28-208-00
Dissection Kit Deluxe
Deluxe Dissecting Kit In A
Zippy Case



S28-209-00
Botany Dissection Kit - Economy
Budget Botany Dissection kit



S28-210-00
DR Botany Dissection Kit
Botany Dissection Kit



S28-211-00
Dissecting Set



S28-212-00
Dissecting Set



S28-213-00
Dissecting Set



S28-214-00
Dissecting Set



S28-215-00
Student Dissecting Set I



S28-216-00
Intermediate Dissecting Set I



S28-217-00
Dissection Premium Kit



S28-218-00
Dissection Premium Kit



S28-231-00
Fish Equram Spatula Set



S28-232-00
Fish Equram Spatula Set



S28-233-00
Fish Equram Spatula Set



S28-251-00
Laboratory Spatula Set

Laborspatel

Espátulas de laboratorio

Spatules de laboratoire

Spatole da laboratorio


S28-252-00

Set of 7 Pcs Spatula stainless steel
 highly flexible polished, wooden handle
 length of blade mm length mm

75	155	16 g
100	190	19 g
130	230	29 g
150	255	34 g
200	320	48 g
250	360	61 g
300	405	61 g


S28-253-00

Set of 5 Pcs
 Spatula stainless steel highly
 flexible plastic handle
 autoclavable
 length of blade mm length mm

75	160	19 g
100	200	26 g
130	240	31 g
150	260	35 g
200	330	63 g


S28-254-00

Spatula stainless steel
 blade flexible
 length of blade 50 mm
 width of blade mm
 15 23 g


S28-255-00

Set of 3 Pcs
 Double Spatula for Mortars
 stainless steel
 width of blade 20 mm rigid
 length mm

150	29 g
210	38 g
250	49 g


S28-256-00

Set of 3 Pcs
 Powder Spatula stainless steel
 width of blade 9 mm
 length mm

150	8 g
185	8 g
210	10 g


S28-257-00

Set of 4 Pcs
 Double Spatula bent
 stainless steel curved
 length mm width of blade mm

150	9	7 g
185	10	9 g
210	12	16 g
250	12	18 g


S28-258-00

Set of 8 Pcs Double Spatulas,



spatula width 15 mm, spatula length 70 mm, length 300 mm
 straight edges, spatula width 6 mm, spatula length 40 mm, length 130 mm
 straight edges, spatula width 6 mm, spatula length 40 mm, length 150 mm
 straight edges, spatula width 6 mm, spatula length 50 mm, length 185 mm
 straight edges, spatula width 6 mm, spatula length 55 mm, length 210 mm
 trapezoid shape, spatula width 5 mm, spatula length 28 mm, length 130 mm
 trapezoid shape, spatula width 5 mm, spatula length 30 mm, length 185 mm
 trapezoid shape, spatula width 5 mm, spatula length 30 mm, length 210 mm


S28-259-00

Set of 6 Pcs
 Micro Double Spatula stainless steel
 flexible
 length mm width of blade mm

100	3	2 g
130	3	2 g
150	3	4 g
185	3	4 g
150	6	5 g
185	6	6 g

Laborspatel, Löffel

Spatules de laboratoire, cuillères

Espátulas de laboratorio, cucharas

Spatole da laboratorio, cucchiari


S28-260-00

Set of 6 Pcs

Double Spatula stainless steel flexible

length mm	width of blade mm
130	9 6 g
150	9 7 g
185	9 8 g
210	11 15 g
250	11 17 g
300	16 42 g


S28-266-00

Micro-spatula set with plastic handle, stainless steel

Set of 4 parts

Content: each 1 spatula 3, 4, 5, 6mm wide Length in mm 150


S28-267-00

Micro-powder Spatula Set with plastic handle, stainless steel

Set of 4 parts Content:

each 1 spatula 3, 4, 5, 6mm wide Length in mm 150


S28-271-00

Set of 6 Pcs Double Spoons,

 Length 150 mm, L x I spoon 1: 11/22 mm, L x I spoon 2: 20/25 mm
 Length 180 mm, L x I spoon 1: 20/27 mm, L x I spoon 2: 24/32 mm
 Length 210 mm, L x I spoon 1: 20/28 mm, L x I spoon 2: 27/35 mm
 Length 250 mm, L x I spoon 1: 24/32 mm, L x I spoon 2: 30/40 mm
 Length 300 mm, L x I spoon 1: 29/40 mm, L x I spoon 2: 35/45 mm
 Length 160 mm, L x I spoon 1: 22/35 mm, L x I spoon 2: 30/45 mm

S28-272-00

Set of 2 Pcs

Double spoons,

length 160 mm, L x I spoon 1: 22/35 mm, L x I spoon 2: 30/45 mm

length 230 mm, L x I spoon 1: 30/60 mm, L x I spoon 2: 38/75 mm


S28-273-00

Set of 2 Pcs

Excavator with button

From Hammacher. Stainless steel Remanit 4024. Autoclavable.

Shovel with button. Length 190 mm. Excavator dimensions: 35/70

Shovel with button. Length 205 mm. Excavator dimensions: 30/63


S28-274-00

Set of 3 Pcs

Micro Spatula Spoon

stainless steel rigid

width of blade 5 mm

width of spoon 5 mm

length mm

150	6 g
185	7 g
210	8 g


S28-275-00

Set of 4 Pcs

Spatula Spoon stainless steel

bent spoon width of spoon 19 mm

length mm

150	19 g
185	23 g
210	28 g
250	33 g

Laborlöffel

Cucharas de laboratorio

Cuillères de laboratoire

Cucchiari da laboratorio


S28-276-00

Set of 4 Pcs

Spatula Spoon stainless steel

rigid broad

length mm width of blade mm

120 23 15 g

150 28 23 g

180 30 35 g

210 30 42 g


S28-277-00

Set of 3 Pcs

Spatula Spoon stainless steel

width of spoon 29 mm

length mm

150 20 g

180 30 g

210 37 g


S28-278-00

Set of 6 Pcs

Double Spoon stainless steel

rigid round handle

length mm width of spoon mm

120 17 / 25 16 g

150 17 / 25 19 g

180 17 / 25 29 g

210 22 / 29 34 g

250 29 / 35 48 g

300 29 / 35 61 g


S28-279-00

Chemist's Spoon

stainless steel

length 150 mm

width 16 mm

12 g


S28-280-00

Spoon, open

stainless steel

length 200 mm

width 25 mm

29 g


S28-281-00

Deflagrating Spoon

cup stainless steel

cup Ø 16 mm

length 450 mm

14 g


S28-282-00

Set of 2 Pcs

Sodium Spoon

copper basket

basket Ø 25 mm

length 180 mm

Single 16 g

double 26 g


S28-283-00

Ladle


stainless steel

cup Ø 75 mm

length 275 mm

73 g

**VATS / MIS THORACOSCOPIC /
MINIMALLY INVASIVE INSTRUMENTS**




LEADER
SURGICAL

Thoracic Scissors




S29-001 Thoracic double jointed curved scissors jaw length 4 cm Working 15,5 cm total length cm 28,5 cm

Thoracic Needle Holder



S29-002 Thoracic Curved Needle Holder jaw length 1,8 cm Working 20 cm total length cm 33 cm

Conventional Scissors




S29-003

MIS Scissors




S29-004

D'Amico Biopsy Forceps



S29-005 6 mm width oval ring jaws Curved width ratchet Working 23 cm total length 33 cm


S29-006 6 mm width oval ring jaws Curved without ratchet Working 23 cm total length 33 cm



S29-007 7 mm width oval cupped jaws Curved width ratchet Working 23 cm total length 33 cm


S29-008 7 mm width oval cupped jaws Curved without ratchet Working 23 cm total length 33 cm

D'Amico Mediastinoscopy Forceps



S29-009 6mm width oval ring jaws curved with ratchet working length 28 cm total length 38 cm


S29-010 6mm width oval ring jaws curved without ratchet working length 28 cm total length 38 cm



S29-011 7 mm width oval cupped jaws curved with ratchet working length 28 cm total length 38 cm

S29-012 7 mm width oval cupped jaws curved without ratchet working length 28 cm total length 38 cm

NEW Rocco Nodule Clamps




S29-013 Curved left ring jaws 3,0 cm inner diameter and 3,7 cm outer diameter with ratchet 24 cm total length 33,5 cm

S29-014 Curved left ring jaws 4,5 cm inner diameter and 5,2 cm outer diameter with ratchet 24 cm total length 33,5 cm

NEW Gonzalez-Rivas Lung Grasping Clamp




NEW Gonzalez-Rivas Lung Grasping Clamp



S29-015 High tension Curved left 20 mm oval ring jaws With ratchet Working 14,5 cm total lenght 24,5 cm

Foerster Lung Grasping Clamps



NEW Foerster Lung Grasping Clamp

S29-016 Curved left 20 mm oval ring jaws With ratchet Working 23,5 cm total lenght 33,5 cm

Foerster Lung Grasping Clamps



Foerster Lung Grasping Clamp


S29-017 Curved right jaws (left-hand style) 12 mm oval ring jaws With ratchet Working 19 cm total length 29 cm

S29-018 Curved left jaws 12 mm oval ring jaws With ratchet Working 19 cm total length 29 cm



Foerster Lung Grasping Clamp

S29-019 Curved left jaws 12 mm oval ring jaws With ratchet Working 23,5 cm total length 33,5 cm



Foerster Lung Grasping Clamps




Foerster Lung Grasping Clamps

S29-020 - Curved right jaws (left-hand style) 20 mm oval ring jaws with ratchet Working 19 cm total length 29 cm

S29-021 - Curved left jaws 20 mm oval ring jaws with ratchet Working 19 cm total length 29 cm

Foerster Lung Grasping Clamps



- S29-022 - Curved left jaws 20 mm oval ring jaws with ratchet Working 24 cm total length 33,5 cm
- S29-023 - Deep Curved left jaws 20 mm oval ring jaws with ratchet Working 24 cm total length 33 cm
- S29-024 - Curved left jaws 25x35 mm oval ring jaws with ratchet Working 24 cm total length 33,5 cm



S29-022




S29-023



S29-024


Duval Lung Grasping Clamps



S29-025 - Curved left jaws 1 cm width triangular ring jaws With ratchet Working 23,5 cm total lenght 34 cm

S29-026 - Curved left jaws 2,5 cm width triangular ring jaws With ratchet Working 23,5 cm total lenght 34 cm

Node Grasping Clamps



Decortication Grasper




S29-030 - Enhanced serrations Curved left jaws With ratchet Working 24 cm total length 33,5 cm

Rivas Dissector



Rivas Dissector




- S29-031 - Mixer-style 1x2 DeBakey jaws With ratchet Working 19 cm total length 29 cm
- S29-032 - Mixer-style 1x2 DeBakey jaws Without ratchet Working 19 cm total length 29 cm
- S29-033 - Mixer-style 1x2 DeBakey jaws Fine tips With ratchet Working 19 cm total length 29 cm
- S29-034 - Mixer-style 1x2 DeBakey jaws Fine tips Without ratchet Working 19 cm total length 29 cm

S29-031 · S29-032 S29-033 · S29-034




DeBakey Cooley Forceps



- S29-035 - DeBakey Cooley jaws With ratchet Working 24 cm total lenght 34 cm
- S29-036 - DeBakey Cooley jaws Without ratchet Working 24 cm total lenght 34 cm

Harken Clamps



- S29-037 - Harken 7.5 cm jaws Tapered 1x2 DeBakey jaws Without ratchet Working 19 cm total length 29 cm
- S29-038 - Harken 7.5 cm jaws Tapered 1x2 DeBakey jaws With ratchet Working 19 cm total length 29 cm
- S29-039 - Harken 7.5 cm jaws Tapered 1x2 DeBakey jaws Without ratchet Working 19 cm total length 29 cm
- S29-040 - Harken 7.5 cm jaws Tapered 1x2 DeBakey jaws With ratchet Working 19 cm total length 29 cm

Harken Clamps




S29-041 - S29-042



S29-043 - S29-044

Harken Clamps




S29-041 - Harken 9.5 cm jaws Tapered 1x2 DeBakey jaws Without ratchet Working 24,8 cm total length 35,5 cm

S29-042 - Harken 9.5 cm jaws Tapered 1x2 DeBakey jaws With ratchet Working 24,8 cm total length 35,5 cm

S29-043 - Harken 9.5 cm jaws Tapered 1x2 DeBakey jaws Without ratchet Working 24,8 cm total length 35,5 cm

S29-044 - Harken 9.5 cm jaws Tapered 1x2 DeBakey jaws With ratchet Working Working 24,8 cm total length 35,5 cm


Dennis Dissector



Dennis Dissector

- S29-045 - Curved left jaws serrated jaws With ratchet Working 23 cm total lenght 33,5 cm
- S29-046 - Curved left jaws serrated jaws Without ratchet Working 23 cm total lenght 33,5 cm


Allis Forceps (Clamp)



Allis Forceps (Clamp)

S29-047 - Curved left jaws 6 mm width jaws With ratchet Working 23 cm total lenght 33,5 cm

D'Amico DeBakey Clamps



S29-048 - Curved left jaws Tapered 1x2 DeBakey jaws 11 cm jaw length Without ratchet Working 25,4 cm total length 37 cm


S29-049 - Curved left jaws Tapered 1x2 DeBakey jaws 11 cm jaw length With ratchet Working 25,4 cm total length 37 cm

S29-050 - Straight jaws Tapered 1x2 DeBakey jaws 11 cm jaw length With ratchet Working 27 cm total length 37,5 cm



S29-049


Chitwood DeBakey Clamps




S29-051 - 6 cm jaw length Curved left jaws 1x2 DeBakey jaws With ratchet Working 21 cm total length 31 cm

S29-052 - 6 cm jaw length Curved left jaws 1x2 DeBakey jaws Without ratchet Working 21 cm total length 31 cm


Chitwood DeBakey Clamps



- S29-053 - 6,5 cm jaw length 1x2 DeBakey jaws Curved left jaws With ratchet Working 23 cm total length 32,5 cm
- S29-054 - 10 cm jaw length 1x2 DeBakey jaws Curved left jaws With ratchet Working 25,4 cm total length 36 cm
- S29-055 - 10 cm jaw length 2x3 DeBakey jaws Curved left jaws With ratchet Working 25,4 cm total length 36 cm



Chitwood DeBakey Clamps



- S29-056 - 7 cm jaw length Curved left jaw 1x2 DeBakey jaws With ratchet Working 17 cm total length 26,5 cm
S29-057 - 7 cm jaw length Curved right jaw 1x2 DeBakey jaws With ratchet Working 17 cm total length 26,5 cm

Super Cut Scissors



S29-058



S29-059



S29-060



S29-061



S29-062



Chitwood Dissecting Scissors

S29-058 - Straight fine blades 1 cm blade length blunt micro tips Working 16,5 cm total length 27 cm

S29-059 - Curved fine blades 1 cm blade length blunt micro tips Working 16,5 cm total length 27 cm



S29-060 - Chitwood dissecting scissors Curved blades 3 cm blade length blunt tips Working 18,5 cm total length 29 cm



Metzenbaum Scissors

S29-061 - Curved blades 3,5 cm blade length blunt tips Working 19,5 cm total length 29 cm

S29-062 - Curved wide blades 3,5 cm blade length blunt tips Working 19,5 cm total length 29 cm



Ring Handle Needle Holders



S29-063



S29-064



S29-065




S29-066




S29-067

Ring Handle Needle Holders



- S29-063 MVT 100.72 - Straight jaws Carbon-inlay jaws Gold ring handle With ratchet Working 16 cm total length 28,5 cm
- S29-064 MVT 100.73 - Curved jaws Carbon-inlay jaws Gold ring handle With ratchet Working 16 cm total length 28,5 cm
- S29-065 MVT 100.74 - Straight jaws Carbon-inlay jaws Gold ring handle With ratchet Working 16 cm total length 28,5 cm
- S29-066 MVT 100.75 - Curved jaws Carbon-inlay jaws Gold ring handle With ratchet Working 16 cm total length 28,5 cm

Ring Handle Needle Holders



- S29-067 - Straight jaws Fine DIAMOND DUST jaws With ratchet Working 16,5 cm total length 27 cm

Ring Handle Needle Holders



S29-068



S29-069



Ring Handle Needle Holders

- S29-068 - Straigth jaws Carbon-inlay jaws Gold ring handle With ratchet Working 21 cm total lenght 34 cm
- S29-069 - Curved jaws Carbon-inlay jaws Gold ring handle With ratchet Working 21 cm total lenght 34 cm



S29-070



S29-071



Ring Handle Needle Holders


- S29-070 - Straigth jaws Carbon-inlay jaws Gold ring handle With ratchet Working 21,5 cm total lenght 34 cm
- S29-071 - Curved jaws Carbon-inlay jaws Gold ring handle With ratchet Working 21,5 cm total lenght 34 cm




Wolf Suction Instruments




Enlarged to
show detail



S29-072



S29-074



S29-077

Variable Suction Regulator
S29-073 - S29-075



- S29-072 - Pediatric Curved blunt tip 6 mm width tip Working 17 cm total length 31 cm
- S29-073 - Pediatric Variable suction Curved blunt tip 6 mm width tip Working 17 cm total length 31 cm
- S29-074 - Curved blunt tip 6 mm width tip Working 22 cm total length 36 cm
- S29-075 - Variable suction Curved blunt tip 6 mm width tip Working 22 cm total length 36 cm
- S29-076 - Curved blunt tip 6 mm width tip Working 22 cm total length 36 cm
- S29-077 - Variable suction Curved blunt tip 6 mm width tip Working 22 cm total length 36 cm

D' Amico Suction Instruments



Enlarged to
show detail

S29-078 - Straight blunt tip 5 mm width tip Working 32 cm total length 45,5 cm

S29-079 - Variable suction Straight blunt tip 5 mm width tip Working 32 cm total length 45,5 cm



Enlarged to
show detail


S29-080 - Straight blunt tip 10 mm width tip Working 32 cm total length 45,5 cm

S29-081 - Variable suction Straight blunt tip 10 mm width tip Working 32 cm total length 45,5 cm


Variable Suction Regulator
S29-080 - S29-081




Dennis Rib Cutter




S29-082 - 10 mm blade 16 mm cutting width Working 29 cm total length 34 cm



Thoracic Tissue Forceps




S29-083 - Thoracic tissue forceps ring jaw length 5 cm Working 18 cm total length cm 33 cm




S29-084 - Thoracic tissue forceps ring jaw length 6 cm Working 18,5 cm total length cm 35,5 cm

Thoracic Vascular Forceps



S29-085 - Thoracic Vascular forceps jaws length 5 cm 30° Working 18 cm total length 33,5 cm



S29-086 - Thoracic Vascular forceps jaws length 6,5 cm 90° Working 18 cm total length 31 cm



LEADER SURGICAL

Manufacturers & Exporters of Surgical

42 A-1 Gulberg 3 Lahore Pakistan

Cell: + 92-322-3000333

Email: info@leadersurgical.org

Website: www.leadersurgical.org

